

D200

STANDARD MANUAL WHEELCHAIR

C86
carbon
grey

- Available in 8 seat widths (from 38 to 52 cm)
- Height and depth adjustable armrests
- Height adjustable push handles
- Inwards and outwards swing-away and removable footrests
- Seat height adjustable in three positions
- Padded seat and backrest
- In height and depth adjustable armpads
- Angle adjustable footrests

The D200 is a lightweight manual wheelchair made of aluminium. Foldable for easy transport and storage without compromising on adjustment possibilities. The new design and carbon colour give it a modern and agreeable look for all.

Features



Height adjustable push handles



Padded seat and backrest



Adjustable front wheels



Two tipping pieces



Easy to fold and compact for transport and storage



Technical information

| | |
|------------------|--|
| ■ Total width | 580 / 600 / 620 / 640 / 660 / 680 / 700 / 720 mm |
| ■ Total length | 1010 mm |
| ■ Total height | 880 - 940 mm |
| ■ Folded width | 300 mm |
| ■ Seat width | 380 / 400 / 420 / 440 / 460 / 480 / 500 / 520 mm |
| ■ Armrest height | 220 - 240 mm |

| | |
|--------------------------|--------------|
| ■ Seat height | 470 - 530 mm |
| ■ Backrest height | 420 mm |
| ■ Seat depth | 430 mm |
| ■ Weight | 16,8 kg |
| ■ Maximum user weight | 130 kg |
| ■ Additional information | |

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.