

Gemini 2

BUGGY



- Head support
- Trunk supports
- Cosy seat with soft fabric
- Attractive design
- Feet stabilization



SCAN the QR code and discover the Gemini 2.

A walk in the park or a family dinner at home, Gemini 2 provides a fully customizable and adjustable solution for the family. Adjusting the seating positioning, changing the color or even adding on an indoor-frame. The Gemini 2 stroller makes it feel like anything is possible!

Thanks to our Easy & Smart positioning control, you can create the correct posture for a child. Even at a young age, perfect seating is important. Therefore Gemini 2 allows for

head, trunk, hip and feet corrections. All made possible by our wide range of positioning accessories and options for this stroller.

2 chairs for 1: Interchange the seat system for a hi-low base frame, changing the Gemini 2 from a stroller into an indoor inclinable chair. Fun for activities at the dinner table and quickly suited for an outdoor trip!



Features



Adjustable pushbar



Footpedal brake

Technical information

■ Total width	620 mm & 690 mm
■ Total height	1200 mm & 1275 mm
■ Total length	950 mm & 1060 mm
■ Folded width	620 mm & 690 mm
■ Seat width	320 mm & 400 mm
■ Height armrests	200 mm

■ Seat height	510 mm
■ Height backrest	520 - 640 mm & 650 - 770 mm
■ Seat depth	300 - 360 mm & 360 - 420 mm
■ Weight	18 kg & 27 kg
■ Maximum user weight	35 kg & 50 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.